

Hybrid-Power

Quick Glance User Guide

drive

DeVilbiss
HEALTHCARE



Need Assistance?
Call 0845 0600 333

1. Patient Weight Settings - 3 comfort levels are available for a variety of patient weights. Use the - / + to adjust the weight settings accordingly.



2. Power - The main power switch is located on the side of the pump. The standby button is located on the front panel.



3. Static Mode - All cells in the mattress will be inflated by pressing the button. To discontinue the action press the alternating button. The mattress will return to the alternating mode after 20 minutes in static.



4. Alternating Mode - One in two cycle. This function will provide optimum pressure relief.



5. Lock - After 5 minutes of inactivity the pump will automatically lock its functions. To unlock the functions press and hold the lock key for 3 seconds.



6. Mute - Press the mute button to silence the alarms.



7. Power Failure / Low Pressure Indicators - When either a power cut or low pressure is detected the pump alarm and respective LED will illuminate.



Hybrid-Power Weight Settings & User Guidelines



Quick Set-up

1. Remove all the existing surfaces from the bed frame.
2. Ensure the printed surface is placed facing upwards.
3. The printed feet should be at the foot end of the bed.
4. Place a loosely fitted sheet over the mattress.
5. Hang the pump at the foot end of the bed; adjust the hangers to secure the pump against the foot board. Avoid placing the pump on the floor.
6. Join the connectors feed pipes to the side of the pump.
7. Plug the pump into the mains electrical outlet.
8. Turn on the pump using the main power switch located on its side. The "Power" button on the pumps front turns off standby mode.

9. The pumps low pressure indicator will remain on until mattress is inflated.
10. Adjust pressure settings according to the patients weight.
11. Once the mattress has inflated, place the patient on the mattress.

Three comfort levels are available for a variety of patient weights. Use the - / + to adjust the weight settings accordingly.

**0-85kg approx
86-170kg approx
171-254kg approx**



**HYBRID POWER
INFORMATION**

Need Assistance?
Call 0845 0600 333

